

*ANNOUNCING the **new PEN** type refractometer*

## **PEN-PRO**



**Price: \$415**

**NEW**

## **PEN-PRO**

**Touch! Dip! One-handed operation!**

**Features - New Continuous Measurement while mixing**

Measurement range	Brix 0.0 to 85.0%
Resolution	Brix 0.1%
Measurement accuracy	Brix +/- 0.2%
Measurement temperature	10 to 60C° (ATC)
Liquid Sample Temperature	10 to 60C°
Ambient temperature	10 to 40C°
Power supply	1 x Size AAA alkaline battery
International Protection class	IP65 (Water Resistant)
Dimensions and weight	16 x 3.8 x 1.8 cm, 70g

- Ø **Water Resistant (IP65)**
- Ø **External Light Interference (ELI)**
- Ø **Automatic Temperature Compensation (ATC)**
- Ø **Calibration with water only**
- Ø **Light & Compact, 70g**

The **PEN-Pro Refractometer** is **EXTREMELY EASY TO USE**

**Digital Hand-held Refractometer . . . from ATAGO-USA**



## Unitech Scientific LLC

### *manufactures wine-testing reagents*

**UniTAB™** and **UniFLEX™**

kits for assessing sugars, organic acids, phenolics and nutrient status

Multi-level Standards

Enzymes & Nucleotide components

### **& also provides**

- ✓ **ChemWell** automated wine analyzer  
& the StatFax analyzers for Enzymatic and ELISA testing
- ✓ **Laboratory & field instrumentation (YSI, Hanna, Atago-USA):**
  - pH Meters
  - Dissolved Oxygen Meters
  - Refractometers
  - Turbidity meters
  - Titration systems
- ✓ **Microbiology media and supplies**
  - **DBDM**, the highly differential Brettanomyces Agar & Broth
  - **Compact Dry**, an innovative dry micro (differential plating) media
  - **Sartorius-Stedim** line of laboratory microbiology and filtration products
- ✓ **Rapid Test devices for Wine**
  - **Z-Brett™**, the 3-hour Brettanomyces kit
  - **HistaSure** 5-minute dip-stick test for Histamine
  - **Quick Test** visually read Wine Test sticks